


- 46600-00**  
Dermatome Handpiece
- Autoclavable
  - 16.000 cpm
  - 300 x 90 x 42 mm / 1170 g



- 46604-00**  
Spare blades  
One side cutting only.  
stainless steel
- 10 ea.
  - 81 x 18,5 x 0,40 mm



- 46603-00**  
Charger unit for 2 battery  
(without battery)
- Mains voltage 100-240 Volts, 50-60 Hz, 0.9 A
  - Charging current 1.5 A
  - Charging voltage 9.6 Volts
  - Protections: Short circuit / Overload / Over voltage
  - Dimensions: 200 x 200 x 119 mm
  - Weight: 2850 g



- 46601-00**  
Rechargeable Ni-MH battery
- Not autoclavable
  - Battery capacity: 1.6 Ah
  - Voltage: 9.6 Volts
  - 136 x 32 x 32 mm / 270 g



- 46602-00**  
Autoclavable sterile funnel  
for aseptic battery loading
- 55 x 55 x 55 mm / 60 g



**46610-00**  
Skin Graft Mesher  
• 205 x 165 x 120 mm / 4000 g



**46611-00** Skin graft carrier > 1,5:1 (Box of 20)  
**46612-00** Skin graft carrier > 3:1 (Box of 20)